EN

1. Intended Use

Disposable Biopsy Forceps are used for taking out biopsy specimens (tissues samples) from the human body, through endoscopes.

2. Product Description

The device is composed by a jaw (in different shapes and sizes), a spiral spring hose (in different lengths and diameters) and an operating handle.

3. Advices before first application

Please read carefully and follow all safety operating instructions and warnings before first application of the device.

A previous knowledge regarding handling and operation is required and essential.

Unpack the instrument carefully and examine the devices for any possibility of damage. In case of any damage or missing items contact your distributor immediately.

If the package present any damaged, the sterility of the device is not guarantee.

4. Handling and Operation

The biopsy forceps must be passed through the endoscope operating channel with its cups closed. Please observe compatibility between the biopsy forceps and the endoscope operating channel diameter.

5. Procedure

- Approach the biopsy forceps near to the tissue (lesion) to be taken from the patient's body.
- Open the cups by pushing the plastic spool at the handle.
- Place the forceps on the tissue and take the sample gently by pulling back the plastic spool (this maneuver will close the forceps cups).
- Keep the forceps cups closed and remove biopsy forceps out of the endoscope carefully.
- This tissue specimen taken will be ready for further histological examination.

6. Storage

This device is delivered sterilized and in order to maintain this state until its first use it must be kept in its original sterilized package.

Do not put any objects on the instrument or its package!

Do not store the instruments near aggressive chemical products!

Do not expose the instruments to direct or indirect sunlight or other ultra-violet rays!

Keep in dry area

Complaints will not be taken under consideration if the instruments have been stored improperly.

7. Cleaning and Sterilization

This device is designed and warranted for single use!

G-Flex Europe S.P.R.L. | 20, Rue de l'industrie 1400 Nivelles - Belgium

Phone: +32(0)67.88.36.65 | Fax: +32(0)67.88.36.88 | Email: sales@g-flex.com | Website: www.g-flex.com

US Commercial Office Only: G-Flex USA, Inc., 74 Southwood Pkwy Suite 600 Atlanta, GA 30354 – USA

EN

8. Risk in case of reuse

This device is for single use only. G-Flex denied all responsibility in case of reuse. The reuse of single use device might present risk to the safety of the user or patient, due to possible uncontrolled contamination and / or unreliability of the device performance.

9. User

The users of G-FLEX instruments must be specialists in their fields. An appropriate and specific training for preparation, care and maintenance of the flexible instruments is required.

10. Repairing and complains

In health protection of our employees only disinfected or sterilized instruments will be accepted for analysis or repairing. In any case disinfection or sterilization dates and validity must be labelled outside the package. If this requirement has not been fulfilled the instruments will be returned to the sender without analysis or repairing.

11. Legal foundation

The Law of the European Union applied.

12. Product support

In case of questions or difficulties concerning our instruments please contact your local distributor or G-FLEX directly during regular working hours.

Times available: Monday to Friday; 9am to 5pm (MET)

13. Symbols

	Date of manufacture		Do not use if package is damaged		Caution
	Use-by date		Keep away from sunlight		Do not reuse
	Batch or lot code		Keep dry		Consult the instruction for use
	Sterilized using ethylene oxide		Temperature limit		Quantity of units per box

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US Commercial Office Only: G-Flex USA, Inc., 74 Southwood Pkwy Suite 600 Atlanta, GA 30354 - USA